



TruCurve Curve Tracker

A First! Track Spinal Curves and Spinal Posture in 20 seconds.




Easy to Operate

1. Turn on device
2. Place at L5
3. Press START
4. Roll to C1
5. Press STOP
6. View results in Patient Report

TruCurve Benefits

- Noninvasive with No Radiation
- Measures Cervical, Thoracic, and Lumbar Curves
- Measures Head Forward and Tilt Forward
- Measures Lateral Tilt
- Displays Saggital and Lateral Views
- TruCurve SmartWord Report captures Six Patient Evaluations
- Ideal for Patient Screenings
- Wireless using standard Bluetooth®
- USB Rechargeable
- Made in the U.S.A.

Simple to read graphs for Saggital and Lateral views. Easily visualize Head Forward and Posture Tilt on the Spinal Graph.



biotools



Made to Measure |

888-850-0550 x111
Olathe, KS 66062
www.biotoolscorp.com